

CP-BP-SFAC

Chemical Properties

CAS No. :

Formula: C₅₀H₃₂N₂O

Molecular Weight:

Keep away from direct sunlight

Storage: Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	CP-BP-SFAC is a luminescent molecule that exhibits intense sky-blue delayed fluorescence when in a pure thin-film state.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481