

TCP-BP-SFAC

## Chemical Properties

CAS No. :

Formula: C74H46N4O

Molecular Weight:

Keep away from direct sunlight

Storage: Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	TCP-BP-SFAC is a luminescent molecule exhibiting intense sky-blue delayed fluorescence in pure thin films. Its photoluminescence (PL) peak is around ~483 nm, with a delayed fluorescence lifetime ranging from 5.4 to 5.7 $\mu$ s.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481