

FITC-C6-LEHD-FMK

Chemical Properties

CAS No. :

Formula: C49H55FN8O14S

Molecular Weight:

Keep away from direct sunlight

Storage: Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	FITC-C6-LEHD-FMK is a fluorescently labeled caspase-9 inhibitor utilized for detecting active caspase-9 in apoptotic mammalian cells. This compound offers a convenient method for the sensitive detection of activated caspase-9 in live cells. Meanwhile, Z-LEHD-FMK serves as a specific caspase-9 inhibitor.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481