


Fluorescent Red Mega 485 NHS-ester

Chemical Properties

CAS No. :	890317-36-9
Formula:	C ₂₉ H ₃₃ N ₃ O ₉ S
Molecular Weight:	599.65
Storage:	Keep away from direct sunlight Store at -20°C <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Fluorescent Red Mega 485 NHS-ester is an amine-reactive fluorescent biolabel suitable for protein labeling.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6676 mL	8.3382 mL	16.6764 mL
5 mM	0.3335 mL	1.6676 mL	3.3353 mL
10 mM	0.1668 mL	0.8338 mL	1.6676 mL
50 mM	0.0334 mL	0.1668 mL	0.3335 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481