

Sulfo-Cyanine5.5 alkyne

Chemical Properties

CAS No. :

Formula: C₄₂H₄₁K₃N₄O₁₃S₄

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Sulfo-Cyanine5.5 alkyne is a near-infrared (NIR) fluorescent dye that features an alkyne-modified Sulfo-Cyanine5.5 (Ex=678 nm, Em=694 nm). It can react with azides for click chemistry labeling.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481