

## ATTO 550 alkyne

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight

Store at -20°C

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	ATTO 550 alkyne is a novel fluorescent label based on the Rhodamine structure, characterized by strong absorption, high fluorescence quantum yield, excellent thermal stability, and photochemical stability. It is suitable for single-molecule detection and high-resolution microscopy. As an alkyne derivative of ATTO 550, it can be used to label proteins or antibodies.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481