

ATTO 532 cadaverin

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**  
**Store at -20°C**  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	ATTO 532 cadaverin is a cadaverine derivative of ATTO 532, used in reactions involving carboxylic acids, aldehydes, and ketones, with maximum excitation/emission wavelengths of 532/552 nm.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481