

ATTO 700

## Chemical Properties

CAS No. :

Formula: C<sub>30</sub>H<sub>36</sub>ClN<sub>3</sub>O<sub>6</sub>S

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	ATTO 700 is an innovative fluorescent dye characterized by a high fluorescence quantum yield, suitable for most labeling applications. It features maximum excitation/emission wavelengths of 700/716 nm.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481