

ATTO 700 NHS ester

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight  
Store at -20°C  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	ATTO 700 NHS ester is a novel fluorescent label based on a Rhodamine structure. It features strong absorption, high fluorescence quantum yield, and excellent thermal and photochemical stability, making it suitable for single-molecule detection and high-resolution microscopy. ATTO 700 NHS ester is a derivative of ATTO 700, capable of labeling proteins or antibodies.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481