

TAMRA-PEG3-NH2

Chemical Properties

CAS No. :

Formula: C34H44N4O7

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	TAMRA-PEG3-NH2 is a TAMRA dye derivative containing three PEG units. It includes an NH2 group, which can undergo a condensation reaction with carboxyl groups to form covalent bonds.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481