

TAMRA-PEG3-Maleimide

Chemical Properties

CAS No. :

Formula: C40H45N5O10

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	TAMRA-PEG3-Maleimide is a TAMRA dye derivative featuring three PEG units. It includes maleimide groups, which are capable of reacting with thiol groups to establish covalent bonds.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481