

## TAMRA-PEG4-NHS

## Chemical Properties

CAS No. :

Formula: C40H46N4O12

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	TAMRA-PEG4-NHS is a derivative of TAMRA dye featuring four PEG units. It contains an NHS group, which can react with amine groups to form an ester bond.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481