

TAMRA-PEG8-NHS

Chemical Properties

CAS No. :

Formula: C48H62N4O16

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	TAMRA-PEG8-NHS is a derivative of TAMRA dye, featuring eight PEG (polyethylene glycol) units. It includes an NHS (N-hydroxysuccinimide) ester group, which can react with amines to form a stable ester bond.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481