

TAMRA-PEG3-COOH

Chemical Properties

CAS No. :

Formula: C34H39N3O9

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	TAMRA-PEG3-COOH is a TAMRA dye derivative that includes three PEG units. Equipped with a carboxyl group, TAMRA-PEG3-COOH can undergo condensation reactions with amines to form covalent bonds.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481