

TAMRA-PEG3-Alkyne

Chemical Properties

CAS No. :

Formula: C₃₄H₃₇N₃O₇

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	TAMRA-PEG3-Alkyne is a TAMRA dye derivative containing three PEG units. It features an alkyne group that can participate in copper-catalyzed azide-alkyne cycloaddition (CuAAC) reactions with molecules possessing an azide group.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481