

## TAMRA-PEG8-Alkyne

## Chemical Properties

CAS No. :

Formula: C44H57N3O12

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	TAMRA-PEG8-Alkyne is a derivative of TAMRA dye with 8 PEG units. It features an Alkyne group that can undergo a copper-catalyzed azide-alkyne cycloaddition reaction (CuAAC) with molecules possessing an Azide group.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481