

ATTO 700 Alkin

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Store at -20°C
Actual storage temperature shall be subject to the COA.

Biological Description

Description	ATTO 700 Alkin is an alkyne derivative of ATTO 700, with maximum excitation/emission wavelengths of 700/716 nm. As a click chemistry reagent, ATTO 700 Alkin contains an Alkyne group that can participate in copper-catalyzed Azide-Alkyne Cycloaddition (CuAAC) reactions with azide-containing molecules.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481