

ATTO 700 Cadaverin

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight
Store at -20°C
Actual storage temperature shall be subject to the COA.

Biological Description

Description	ATTO 700 Cadaverin is a derivative of ATTO 700 and is utilized in carboxylic acid, aldehyde, and ketone reactions, with maximum excitation/emission wavelengths of 700/716 nm.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481