

ATTO 647 maleimide

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Store at -20°C
Actual storage temperature shall be subject to the COA.

Biological Description

Description	ATTO 647 Maleimid is a maleimide ester derivative of ATTO 647, capable of labeling proteins or antibodies, with maximum excitation/emission wavelengths of 630/651 nm.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481