

## ATTO 565 alkyne

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Keep away from direct sunlight

Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	ATTO 565 alkyne is a novel fluorescent label based on the Rhodamine structure. It possesses strong absorption, high fluorescence quantum yield, and excellent thermal and photochemical stability, making it suitable for single-molecule detection and high-resolution microscopy. As an alkyne derivative of ATTO 565, it can be used to label proteins or antibodies.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481