

## ATTO 488 alkyne

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**  
**Store at -20°C**  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	ATTO 488 alkyne is a novel fluorescent label based on the Rhodamine structure, featuring strong absorption, high fluorescence quantum yield, superior thermal stability, and photochemical stability. It is ideal for single-molecule detection and high-resolution microscopy. As an alkyne derivative of ATTO 488, it can be used to label proteins or antibodies.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

**This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use**

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481