

ATTO 665 maleimid

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Store at -20°C
Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|--|
| Description | ATTO 665 maleimid is a maleimide ester derivative of ATTO 665, used for labeling proteins or antibodies, with maximum excitation/emission wavelengths of 663/680 nm. |
| Targets(IC50) | Others |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481