

FAM-SAMS

Chemical Properties

CAS No. :

Formula: C₉H₁₄N₂O₂S₂

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	FAM-SAMS is a 5-FAM-labeled SAMS peptide, serving as a specific substrate for AMP-activated protein kinase (AMPK).
Targets(IC50)	AMPK

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481