

iso-VQA-ACC TFA

## Chemical Properties

CAS No. :

Formula: C<sub>29</sub>H<sub>35</sub>N<sub>7</sub>O<sub>9</sub>.xC<sub>2</sub>H<sub>5</sub>F<sub>3</sub>O<sub>2</sub>

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	iso-VQA-ACC TFA is a substrate of the constitutive proteasome.
Targets(IC50)	Proteasome

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481