

D-Val-Gly-Arg-pNA

Chemical Properties

CAS No. :

Formula: C19H29N7O6

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	D-Val-Gly-Arg-pNA serves as a chromogenic substrate for TPA (tissue plasminogen activator). It can be utilized to assess the amidolytic activity of both TPA I and TPA II.
Targets(IC50)	Serine Protease

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481