

## Penetratin-PI3K $\gamma$ (126-150)

### Chemical Properties

CAS No. :

Formula: C229H362N76O57S

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Penetratin-PI3K $\gamma$ (126-150) is a peptide inhibitor of PI3K $\gamma$ that plays a significant role in respiratory diseases.
Targets(IC50)	PI3K

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481