

APT STAT3, scrambled

Chemical Properties

CAS No. :

Formula: C137H190N42O43

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	APTSTAT3, scrambled (APTscr) is a control peptide for the STAT3-specific peptide (APTSTAT3) and does not bind STAT3. APTSTAT3, scrambled retains the same trpzip scaffold, but features a disordered sequence at the target binding site.
Targets(IC50)	STAT

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481