

retro-inverso TAT-Beclin 1 (D-amino acid)

## Chemical Properties

CAS No. :

Formula: C164H251N57O45

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	The compound retro-inverso TAT-Beclin 1 D-amino acid demonstrates superior activity and resistance to proteolytic degradation compared to the L-amino acid peptide when used in vivo. TAT-Beclin 1 is capable of inducing autophagy in peripheral tissues of adult mice and in the central nervous system of newborn mice.
Targets(IC50)	Autophagy

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481