

APTSTAT3-9R, scrambled

## Chemical Properties

CAS No. :

Formula: C202H315N83O58

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

|               |  |
|---------------|--|
| Description   | APTSTAT3-9R, scrambled (APTscr-9R), is a control peptide structurally similar to APTSTAT3-9R, but features a disordered sequence in the target binding domain. |
| Targets(IC50) | STAT   |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481