

Compstatin control peptide TFA

Chemical Properties

CAS No. :

Formula: C66H101N23O17.xC2HF3O2

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Compstatin control peptide TFA serves as a negative control for Compstatin, a potent inhibitor of the complement system's C3 (complement system C3).
Targets(IC50)	Complement System

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481