

Nuclear pore complex protein Nup98 (315-360)

Chemical Properties

CAS No. :

Formula: C190H288N50O64S

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Nuclear pore complex protein Nup98 (315-360) represents the segment of the nuclear pore complex (NPC) protein ranging from residues 315 to 360.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481