

Ac-VDID TFA

## Chemical Properties

CAS No. :

Formula: C<sub>23</sub>H<sub>35</sub>F<sub>3</sub>N<sub>4</sub>O<sub>12</sub>

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

|               |  |
|---------------|--|
| Description   | Ac-VDID TFA is a short peptide sequence with an acetyl group (Ac) at the terminal end. |
| Targets(IC50) | Others   |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481