

G12 TFA

Chemical Properties

CAS No. :

Formula: C54H96F3N15O17

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	G12 (Ras5-17) TFA is a wild-type Ras peptide comprising amino acids 5-17 (KLVVVGAGGVGKS). It serves as a control in studies involving mutant Ras peptides, such as V12.
Targets(IC50)	Ras

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481