

## Neurogranin (48-76), human

### Chemical Properties

CAS No. :

Formula: C97H161N39O34

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Neurogranin (48-76), human is a predominant endogenous peptide found in the brain tissue of Alzheimer's disease patients. In studies related to Alzheimer's disease (AD), Neurogranin (48-76), human shows potential as a biomarker for synaptic function.
Targets(IC50)	CaMK,PKA,PKC

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481