

Phenylpyruvate tautomerase

Chemical Properties

CAS No. : 9023-54-5

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Phenylpyruvate tautomerase catalyzes the interconversion between the enol and keto isomers of phenylpyruvate and (para-hydroxyphenyl)pyruvate.
Targets(IC50)	Endogenous Metabolite

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481