

Cy3-PEG3-endo-BCN

Chemical Properties

CAS No. :

Formula: C₄₉H₆₇ClN₄O₆

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cy3-PEG3-endo-BCN is a derivative of the dye Cyanine 3 (Cy3) containing three PEG units. It includes a lipophilic bidentate macrocyclic ligand (BCN), which can be utilized for synthesizing macrocyclic complexes. In click chemistry, endo-BCN can form stable triazoles with azide-containing molecules under catalyst-free conditions.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481