

Cy3-PEG2-SCO

Chemical Properties

CAS No. :

Formula: C₄₅H₆₁ClN₄O₅

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cy3-PEG2-SCO is a variant of the Cyanine 3 (Cy3) dye containing two PEG units. It features an SCO group, which can form covalent bonds with amino groups. This SCO group is typically employed to react with amino acid residues in proteins or peptides, particularly lysine.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481