

Cy3-PEG7-SCO

Chemical Properties

CAS No. :

Formula: C₅₅H₈₁ClN₄O₁₀

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cy3-PEG7-SCO is a Cyanine 3 (Cy3) dye derivative containing seven PEG units. This compound features an SCO group, enabling covalent bonding with amino groups. SCO is commonly employed to react with amino acid residues of proteins or polypeptides, particularly lysine.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481