

Cy5-PEG3-Tetrazin

Chemical Properties

CAS No. :

Formula: C₅₁H₆₅ClN₈O₅

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cy5-PEG3-Tetrazin is a derivative of Cyanine 5 (Cy5) dye comprising three PEG units. It features a Tetrazine group that can engage in an inverse electron demand Diels-Alder reaction (IEDDA) with molecules containing a TCO group.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481