

Cy3-PEG3-Alkyne

Chemical Properties

CAS No. :

Formula: C₃₉H₅₂ClN₃O₄

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cy3-PEG3-Alkyne is a derivative of the Cyanine 3 (Cy3) dye that incorporates three PEG units. It features an alkyne group capable of undergoing a copper-catalyzed azide-alkyne cycloaddition (CuAAC) reaction with azide-containing molecules.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481