

Cy3-PEG4-Alkyne

Chemical Properties

CAS No. :

Formula: C₄₁H₅₆ClN₃O₅

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cy3-PEG4-Alkyne is a derivative of the Cyanine 3 (Cy3) dye that includes four PEG units. It features an alkyne group, facilitating copper-catalyzed azide-alkyne cycloaddition (CuAAC) with molecules that possess an azide group.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481