

## Cy5-PEG2-TCO

### Chemical Properties

CAS No. :

Formula: C47H65ClN4O5

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Cy5-PEG2-TCO is a fluorescent dye that reacts with tetrazine-containing molecules. It features a TCO group, enabling it to participate in an inverse electron demand Diels-Alder reaction (IEDDA) with tetrazine groups.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481