

## Cy5-PEG3-TCO

### Chemical Properties

CAS No. :

Formula: C<sub>49</sub>H<sub>69</sub>ClN<sub>4</sub>O<sub>6</sub>

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Cy5-PEG3-TCO is a Cy5 dye derivative containing three PEG units. The TCO group in Cy5-PEG3-TCO enables it to undergo an inverse electron-demand Diels-Alder reaction (IEDDA) with molecules bearing a Tetrazine group.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481