

Cy5-PEG7-TCO

Chemical Properties

CAS No. :

Formula: C₅₇H₈₅ClN₄O₁₀

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cy5-PEG7-TCO is a derivative of Cyanine 5 (Cy5) dye containing seven PEG units. It can undergo an inverse electron-demand Diels-Alder reaction (iEDDA) through its TCO group with molecules that have Tetrazine groups.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481