

Cy5-PEG3-SCO

Chemical Properties

CAS No. :

Formula: C₄₉H₆₇ClN₄O₆

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cy5-PEG3-SCO is a Cyanine 5 (Cy5) dye derivative containing three PEG units. The compound features an SCO group that enables covalent bonding with amino groups. Typically, SCO is utilized for reactions with amino acid residues of proteins or peptides, especially lysine.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481