

IMP-2373

## Chemical Properties

CAS No. :

Formula: C<sub>21</sub>H<sub>20</sub>N<sub>4</sub>O<sub>2</sub>

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	IMP-2373 is the first activity-based probe (ABP). It is a covalently conjugated pan-deubiquitinase (DUB) ABP designed to monitor DUB activity within living cells.
Targets(IC50)	DUB

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481