

AgAuSe-COOH (1000 nm)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Store at -20°C
Actual storage temperature shall be subject to the COA.

Biological Description

Description	AgAuSe-COOH (1000 nm) are fluorescent quantum dots that emit fluorescence in the near-infrared II region, with an emission peak at 1000 nm. These dots possess advantages such as a narrow bandgap, large Stokes shift, and excellent photostability. AgAuSe-COOH (1000 nm) can bind to proteins, antibodies, peptides, and PEG, making them versatile for use in the field of bioimaging.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481