

## Sulfo-Cy5.5 dUTP

### Chemical Properties

CAS No. :

Formula: C58H69N6O28P3S4

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Sulfo-Cy5.5 dUTP is a dUTP derivative with far-red fluorescent properties, featuring a sulfonated Cyanine5.5 group.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481