

Endoproteinase Asp-N

Chemical Properties

CAS No. : 55576-49-3

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Endoproteinase Asp-N (Asp-N) is a metalloprotease that specifically cleaves on the N-terminal side of aspartic acid and sulfoalanine residues.
Targets(IC50)	MMP

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481