


## Goserelin EP Impurity E

## Chemical Properties

CAS No. :	147688-42-4
Formula:	C58H83N17O13
Molecular Weight:	1226.39
Storage:	Keep away from moisture Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	GoserelinEP Impurity E is an impurity found in Goserelin.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8154 mL	4.077 mL	8.154 mL
5 mM	0.1631 mL	0.8154 mL	1.6308 mL
10 mM	0.0815 mL	0.4077 mL	0.8154 mL
50 mM	0.0163 mL	0.0815 mL	0.1631 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481